

# QAZ-1



Sample Splitter  
**QAZ-1**

**ZACCARIA**

# QAZ-1

**Model:**

QAZ-1  
.....

**Sample capacity:**

Product density: 1 kg = 800 kg/m<sup>3</sup> or 0.8 kg/l  
.....

**Approx. shipping volume (dcm<sup>3</sup>):**

60  
.....

**Approx. shipping weight (kg):**


11  
.....

## Sample Splitter QAZ-1

Zaccaria with its complete range of rice equipment, introduce the newest equipment in the laboratory line: Sample Splitter like "Jones".

This equipment is developed to split off, standardization and reduction of cereals samples, granulate product and others.

Developed in stainless steel with twenty canals of 12mm of width, disposed side by side in alternate way, where ensures high precision on the splitting work. Another important characteristic is the construction is developed by boxes and avoid obstruction points for the product pass.



OUR EXPERIENCE MAKES THE DIFFERENCE

**INDÚSTRIAS MACHINA ZACCARIA S/A**

Rua Laranjal, 180  
CEP: 13484-016 - Limeira - SP - Brasil  
Phone: (55.19) 3404 5725  
Fax: (55.19) 3404 5720

[www.zaccaria.com.br](http://www.zaccaria.com.br)  
[comex@zaccaria.com.br](mailto:comex@zaccaria.com.br)



.....  
Indústrias Machina Zaccaria S/A reserves the right to change the information contained in this catalog, including the color and details of the machine, without previous notice.